

## Protein Extraction Buffer N

**Designed for native applications preserving protein conformation and biological activity**

### INTRODUCTION

The Covaris **Protein Extraction Buffer N** (PN 520099) is a physiological buffer optimized for the extraction of proteins in their native conformation. The reagent does not use polyethylene oxide series detergents such as Triton X-100 for protein extraction, since the detergent has a micellar molecular weight of 80,000 Daltons and is difficult to remove by dialysis. Native proteins extracted in the Protein Extraction Buffer N are directly compatible with immunoassay, protein assay, and Blue Native PAGE. Since the reagent does not diminish protein binding, it can be used directly with most affinity enrichment schemes.

### KIT CONTENTS

- Bottle N ..... Protein Extraction Buffer N, 18 mL
- Vial H ..... 100X Protease Inhibitor Cocktail, 0.25 mL (optional)
- EDTA ..... 100X EDTA, 0.25 mL (optional)

### STORAGE

Kit components can be stored at 4°C for one year.

### PREPARATION OF KIT REAGENTS

1. Bottle N is ready to use. Do not dilute prior to use.
2. When applicable, add Protease Inhibitor Cocktail (Vial H) to the prepared Buffer N immediately before use. **NOTE: Some protease inhibitors are irreversible and may have deleterious effects on biological activity. Add protease inhibitors when it is determined they will not interfere with downstream analysis.**
3. When applicable, add EDTA (Vial E) to the prepared Buffer N. **NOTE: EDTA reversibly inactivates enzymes by the chelation of metal ion cofactors. EDTA has been shown to denature some proteins.**

### PROTEIN EXTRACTION PROTOCOL

The below protocol is abbreviated. For detailed protein extraction protocols with cryoPREP Extraction Systems, including the t-PREP, please refer to <http://covarisinc.com/resources/protocols/>.

1. Add the appropriate volume of Native Protein Extraction Buffer N to cells or cryofabricated tissue in the respective AFA Tube or t-PREP for homogenization in a Focused-ultrasonicator.
2. Process samples in a Covaris Focused-ultrasonicator according to the conditions listed in Table 1. A time course should be conducted to determine the optimum conditions for each tissue type and protein of interest.
3. Following AFA, transfer the sample to a 1.5 mL microcentrifuge tube and centrifuge at maximum speed to pellet cellular debris. Transfer the supernatant to a new 1.5 mL microcentrifuge tube for analysis.
4. Buffer N is compatible with the Bradford protein assay and other similar protein assays to determine the protein concentration. Buffer N should be used as the assay blank and for serial dilution of protein standard.

	TC12x12	t-PREP
Volume	1 mL	250 mL
Tissue Mass	< 200 mg	< 10 mg
Duty Cycle	10%	5%
Peak Incident Power (S220 and E220) Intensity (S2 and E210)	70 Watts Level 2	35 Watts Level 1
Cycles per burst	200	
Processing Time	Empirical	
Bath Temperature	7°C	
Power Mode (S2 and E210)	Frequency sweeping	
Degas Mode	Continuous	
Water level (RUN)* S2/S220 E210/E220	Level 8 Level 5	Level 5 NA

